

NEOCERA Chamber Standard Operating Procedure (SOP)

Loading/Unloading

1. Set N₂ pressure to about 10kPa
2. Press ↓ on cold cathode controller to turn off gauge
3. Close gate valve
4. Flip **Pump** switch to turn off RP and TMP
5. Open black valve on N₂ line (near regulator)
 - a. Open process gas valve on back side of the chamber
6. Loosen nuts on port
7. Wait for TMP to reach ~500Hz. Solenoid valve will open (should be able to hear) to vent TMP.
8. Open **Vent** valve
 - a. Note: The pressure needs to be at atmosphere within two minutes.
9. Close the black valve, the **Vent** valve, and the process gas valve once atmosphere is reached

10. Disconnect power and thermocouple cables and remove heater column
 - a. Make sure to cover opening with white plastic cover once heater column is removed
11. Load clean substrate onto sample holder by tightening clips over top the substrate.

Pumping Down

1. Insert heater column and secure with nuts
2. Make sure **Process Gas** and **Course** valves are closed
 - a. Open **Process Gas** only, if you're going to grow in process gas
3. Open gate valve
4. Flip **Pump** switch and wait for pressure to drop to ~1 mTorr
5. Press ↑ on cold cathode gauge controller
6. Wait ~30 minutes to reach 10⁻⁶ Torr

Deposition

1. Move middle mirror in optics tower to zero.
2. Open **Control Target Motors** window in the Neocera PLD program.
3. Use single laser pulses to figure out where to position the turret
4. Use both rotation and rastering (generally target position +- 5-10 degrees)
5. Turn on process gas flow if you're using one.
 - a. Allow to flush for at least 5 minutes
6. Turn on the power to the temperature controller and the heater. Set to desired substrate temperature and set ramp rate in Neocera PLD program.
 - a. Don't exceed 10C/min
 - b. Ramp up and ramp down can be programmed by the computer
8. After finished depositing, turn off raster and rotation.
9. Shut off process gas if needed.

Process Gas

1. Make sure **Process Gas** valve is open before you start pumping down. Gas cylinder should be closed.
2. Pump down the chamber.
3. Turn switches on MFC to **On** and **Flow**
4. Adjust flow rate to get desired pressure
 - a. If going above 10 mTorr, almost close the gate valve.